VERMEIREN

Kim

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI

EN Instructions to the specialist dealer

This instruction manual is part and parcel of the product and must accompany every product sold.

Version: C, 2014-05

All rights reserved, including translation.

No part of this manual may be reproduced in any form what so ever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated or distributed by using electronic systems.

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© N.V. Vermeiren N.V. 2014

EN Contents

1. Technical details	2
2. Components	2
3. Explanation of symbols	2
4. Check after delivery	2
5. Intended use	2
6. Assembly / Dismantling	2
6.1. Assembly pedal exerciser	2
6.2. Dismantling pedal exerciser	2
7. Folding / Unfolding pedal exerciser	3
7.1. Unfolding pedal exerciser	3
7.2. Folding pedal exerciser	3
8. Mounting calorie meter	3
9. Working calorie meter	3
10. Exercises of the legs	4
11. Resistance adjustment	4
12. Safety instructions	4
13. Maintenance, inspection, cleaning	4

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Pedal Exerciser Kim

1. Technical details

⚠ WARNING: Risk of unsafe settings - Use only the settings described in this manual.

Make	Vermeiren
Address	Vermeirenplein 1/15, B-2920 Kalmthout
Type	Pedal exerciser
Model	Kim 
Width	400 mm
Length	500 mm
Height	300 mm
Folded width	400 mm
Folded length	300 mm
Folded height (calorie meter, adjustable star button removed, legs folded)	170 mm
Weight	2,35 kg
Chargeable weight	max. 50 kg
Storage temperature	+5°C - +41°C
We reserve the right to introduce technical changes. Measurement tolerance +/- 15 mm / 1,5 kg	

First of all we wish to thank you for the trust you placed in us by selecting a **VERMEIREN** product.


Before using the pedal exerciser, read the instruction manual carefully: it will familiarise you with this product.

If you have any further questions, please consult your specialist dealer.

To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility. A list of Vermeiren facilities can be found on the last page.

Visually impaired people can contact the dealer for the instructions for use.

2. Components

3. Explanation of symbols

Maximum mass

CE conformity

Type designation

4. Check after delivery

Unpack your product and verify if the delivery is complete. Following items should be included:

- pedal exerciser frame with 2 pedals
- Adjustable star button (resistance)
- Calorie meter (without battery)
- manual

Verify your product for transport damage. If you find any damages after delivery proceed as follows:

- contact the transporter
- have a list made of any problems
- contact your supplier

5. Intended use

The pedal exerciser shall be used by 1 person.

The pedal exerciser shall be used indoor.

The construction of the pedal exerciser allows full use by elderly people or during revalidation in home environment.

The pedal exerciser can increase the flexibility by exercising the range of motion of the joints of the lower body (exercise of the legs).

The pedal exerciser shall not be used as toy for children.

The pedal exerciser should only be used on horizontal, dry and non smooth surfaces.

The manufacturer is not liable for damage caused by the lack of or improper service or as a result of not following instructions from this manual.

6. Assembly / Dismantling

6.1. Assembly pedal exerciser

- Unfold the pedal exerciser.
- Mount the adjustable star button on the pedal exerciser.
- Mount the calorie meter on the pedal exerciser.


6.2. Dismantling pedal exerciser

- Remove the adjustable star button from the pedal exerciser.
- Remove the calorie meter from the pedal exerciser.
- Fold the pedal exerciser.

7. Folding / Unfolding pedal exerciser

- ⚠ **CAUTION:** Risk of clamping - Keep fingers, cloths, jewelry or items away from moving parts of the pedal exerciser (folding legs and upper frame when folding the pedal exerciser).

To transport the pedal exerciser in the car or for storage of the pedal exerciser you can fold the pedal exerciser.


7.1. Unfolding pedal exerciser

- Fold the legs ② outwards until the locking pin ① clicks automatically in the hole.
- Check that the legs ② are firmly fixated and cannot move.

7.2. Folding pedal exerciser


- Press the locking pin ①.
- Fold the legs inwards ②.

8. Mounting calorie meter

- Place the calorie meter ① on the plate ② of the pedal exerciser.

9. Working calorie meter

- Press the red selection button to activate the calorie meter or start pedaling with your pedal exerciser. When the magnet pass by the calorie meter the calorie meter activates automatically.
- The calorie meter display turns off automatically after 4 à 5 minutes.
- When pedaling with the pedal exerciser the STOP in the upper left corner disappears, and all the functions starts with calculating.
- If you stop pedaling after 4 seconds the STOP will reappears in the upper left corner, all functions stop with calculating.
- The calorie meter gives 4 different functions.

CNT: Gives the amount of revolutions of the pedals (count how many the magnet on the pedal pass the calorie meter) from these exercise session.

TOT.CNT: Gives the total amount of revolutions of the pedals from different exercise sessions.

CAL: Gives the total calories amount consumed from the exercise session.

TIME: Gives the total time of the exercise session.

Press the red selection button to go to the next function (TOT.CNT, CAL, TIME). The display show always the CNT function.

Press the red selection button for 4 seconds to reset the functions (CNT, CAL, TIME). The TOT.CNT is not reset. When you have reset the data a next exercise session shall be started.

- **SCAN:** Press the red selection button to go to scan. When the arrow pointing to SCAN you have selected the mode SCAN. Now the display goes automatically through all the functions (TIME, TOT.CNT, CAL), and shows them consecutive on the display.
- When you remove the battery all the data (TIME, TOT.CNT, CNT, CAL) shall be lost. If you see an improper display, please replace the battery to have a good result. Use a 1,5 V battery as power supply.

10. Exercises of the legs

- ⚠ **CAUTION: Risk of injury - Sit always in a chair when performing the leg exercise. Do not stand on the pedal exerciser.**
- ⚠ **CAUTION: Risk of injury - Make sure the chair does not move or rotate. If the chair have wheels, be sure the brakes are engaged.**
- ⚠ **CAUTION: Risk of injury - Be sure you wear good shoes that not slip off the pedals. (Rubber soled shoes)**
- ⚠ **CAUTION: Risk of burns – Do not take the pedal exerciser to the pedals tube after training, this could be warm caused of the friction of moving parts.**
 - Place pedal exerciser on a horizontal, non smooth surface directly in front of the chair.
 - Sit comfortable in a non movable chair.
 - Place your feet in the pedals.
 - Now you can pedal with the pedal exerciser.

11. Resistance adjustment

Turn on the adjustable star button to adjust the resistance of pedaling.

Turn clockwise (direction +):

Resistance increase -- heavier pedaling.

Turn counterclockwise (direction -):

Resistance decrease -- looser pedaling.

12. Safety instructions

- Consult your specialist (doctor) before using the pedal exerciser.
- Stop the exercises and consult your specialist immediately if you feel

dizziness, severe muscle/joint pain or chest pain.

- Warm-up and stretch about 5 minutes before beginning your exercises.
- Exercises shall be in a controlled manner and begin slowly.
- Do not exhaust you.
- **Do not exceed the chargeable weight. The chargeable weight for your pedal exerciser is 50 kg.**
- Use the pedal exerciser on horizontal, dry, non-smooth and indoor surfaces.
- Keep a minimum safety area clearance of 1 meter of the pedal exerciser. Be sure these area remains clear during use of the pedal exerciser.
- Keep it away from children. The pedal exerciser shall not be used as a toy.
- When it is used by, on or near children, invalids or disabled persons close supervision is necessary.
- Risk of burns – Be careful when using in hot or cold environments (sunshine, extreme cold, saunas, etc.) for a sufficient amount of time and when touching.
- Do not place fingers between components of the pedal exerciser, when folding.
- Keep fingers, cloths, jewelry or items away from moving parts of the pedal exerciser (rotating pedals).
- Wear comfortable and suitable clothing when using the pedal exerciser.
- Be sure you wear good shoes that not slip off the pedals.
- Make sure the chair where you sits in does not move or rotate.
- Do not stand on the pedal exerciser.
- Do not use the pedal exerciser when it is not working properly.
- Take note of the instructions for care and service. The manufacturer is not liable for damage caused by improper servicing/care.

13. Maintenance, inspection, cleaning

For the maintenance manual of the pedal exerciser refer to the Vermeiren website: www.vermeiren.be.

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ENGLISH

CONTRACTUAL WARRANTY

We offer 5 years of warranty on standard wheelchairs, lightweight wheelchairs 4 years. Electronic wheelchairs, tricycles, beds and other products: 2 years (batteries 6 months) and multiposition wheelchairs 3 years. This warranty is limited to the replacement of defective or spare parts.

APPLICATION CONDITIONS

In order to claim this warranty, part "B" of this card has to be given to your official Vermeiren dealer. The warranty is only valid when parts are replaced by Vermeiren in Belgium.

EXCEPTIONS

This warranty is not valid in case of:

- damage due to incorrect usage of the wheelchair,
- damage during transport,
- involvement in an accident,
- a dismount, modification or repair carried outside of our company and/or official Vermeiren dealership,
- normal wear of the wheelchair,
- non-return of the warranty card

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N.V. VERMEIREN N.V.
Vermeirenplein 1/15
B-2920 Kalmthout
BELGIUM

Notes

A series of horizontal dotted lines for writing notes.

Notes

A series of horizontal dotted lines for writing notes.

Belgium

N.V. Vermeiren N.V.

Vermeirenplein 1 / 15
B-2920 Kalmthout
Tel: +32(0)3 620 20 20
Fax: +32(0)3 666 48 94
website: www.vermeiren.be
e-mail: info@vermeiren.be

France

Vermeiren France S.A.

Z. I., 5, Rue d'Ennevelin
F-59710 Avelin
Tel: +33(0)3 28 55 07 98
Fax: +33(0)3 20 90 28 89
website: www.vermeiren.fr
e-mail: info@vermeiren.fr

Italy

Reatime S.R.L.

Viale delle Industrie 5
I-20020 Arese MI
Tel: +39 02 99 77 07
Fax: +39 02 93 58 56 17
website: www.reatime.it
e-mail: info@reatime.it

Poland

Vermeiren Polska Sp. z o.o

ul. Łączna 1
PL-55-100 Trzebnica
Tel: +48(0)71 387 42 00
Fax: +48(0)71 387 05 74
website: www.vermeiren.pl
e-mail: info@vermeiren.pl

Czech Republic

Vermeiren ČR S.R.O.

Nadrazni 132
702 00 Ostrava 1
Tel: +420 596 133 923
Fax: +420 596 133 277
website: www.vermeiren.cz
e-mail: info@vermeiren.cz

Germany

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Wahlerstraße 12 a
D-40472 Düsseldorf
Tel: +49(0)211 94 27 90
Fax: +49(0)211 65 36 00
website: www.vermeiren.de
e-mail: info@vermeiren.de

Austria

L. Vermeiren Ges. mbH

Winetzhammerstraße 10
A-4030 Linz
Tel: +43(0)732 37 13 66
Fax: +43(0)732 37 13 69
website: www.vermeiren.at
e-mail: info@vermeiren.at

Switzerland

Vermeiren Suisse S.A.

Hühnerhubelstraße 59
CH-3123 Belp
Tel: +41(0)31 818 40 95
Fax: +41(0)31 818 40 98
website: www.vermeiren.ch
e-mail: info@vermeiren.ch

Spain / Portugal

Vermeiren Iberica, S.L.

Carratera de Cartellà, Km 0,5
Sant Gregori Parc Industrial Edifici A
17150 Sant Gregori (Girona)
Tel: +34 972 42 84 33
Fax: +34 972 40 50 54
website: www.vermeiren.es
e-mail: info@vermeiren.es